USER MANUAL

---

**OAKWORKS®**

**Imaging and Pain Management  
Accessory - Carbon Fiber Armboard**

# USER MANUAL - SINGLE CARBON FIBER ARM BOARD

## MODELS

Standard Arm Board:

- 3395-01: with 1" pad (2.5cm)
- 3395-05: with 2" pad (5cm)
- 3395-07: with 3" pad (7.5cm)
- 3395-06: with 4" pad (10cm)
- 3395-04: with 2" & 4" pads (5cm & 10cm)

## PART IDENTIFICATION

Standard Arm Board Pads:

- 4405-06: 1" pad (2.5cm)
- 4407-06: 2" pad (5cm)
- 4409-06: 3" pad (7.5cm)
- 4408-06: 4" pad (10cm)

## USE

This radiolucent device can be used with any flat support surface where ionizing radiation (x-ray) is utilized for imaging. This arm board can be used for imaging in any environment including private practice, hospital & surgery. It can be used to support arms and feet.

*Note: This product is an imaging platform. It is not a surgical platform*

## DIRECTIONS

Unfold the Arm Board and slide the **Base Section** under the patient and the table top. Position the patient and the **Support Section** to suit the needs of the procedure. The weight of the patient will hold the Arm Board in place.


Rotate the Support Section to the desired position. This section rotates up to 180°.

Prone Position at 135°

Supine Position at 180°

### **CAUTION**

DO NOT exceed 30 lbs. (13.6 kg) of load on the **Support Section**. Injury can occur.

To remove the Arm Board, ask the patient to move their limb off the **Support Section** and lift their body slightly. This will allow you to pull out the **Base Section** of the Arm Board. Fold the **Base and Support Sections** for storage.


## CLEANING & DISINFECTION

Reference the Recommended Cleaners and Disinfectant list (MMINML0008-EN) that came with the product or on our website [www.oakworksmmed.com](http://www.oakworksmmed.com) under product information.

## SPECIFICATIONS

Dimensions (Width x Length) Standard Arm Board	Base Section: 7.5" x 24" (19cm x 61cm) Support Section: 6" x 23" (15cm x 58cm)
Maximum Weight Rating	30 lbs. (13.6 kg)
Radiolucency	Arm Board: Maximum of 1.30mm AL @ 100 kVp. HVL of 3.6mm 1" Pad (2.5cm): Maximum of .42mm AL @ 100 kVp. HVL of 3.6mm 2" Pad (5cm): Maximum of .57mm AL @ 100 kVp. HVL of 3.6mm 3" Pad (7.5cm): Maximum of .68mm AL @ 100 kVp. HVL of 3.6mm 4" Pad (10cm): Maximum of .79mm AL @ 100 kVp. HVL of 3.6mm

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# UNIQUE DEVICE IDENTIFICATION (UDI) INFORMATION

## SYMBOL IDENTIFICATION

This symbol, when used in this manual and on product labels, represents a caution warning. Be sure to read and comply with all precautions and warnings.

This symbol, when used in this manual and on product labels, indicates the name and address of the manufacturer.

This symbol, when used in this manual or on product labels, indicates the country of manufacture along with date of manufacture of the device next to it.

## UNIQUE DEVICE IDENTIFICATION (UDI) INFORMATION

The UDI Label is located on the underside of the arm board

### CARBON FIBER ARM BOARD

**Carbon Fiber Arm Board-2in Pad**

<b>Manufacturer</b>	<b>GTIN</b>	00817463022506
OAKWORKS, Inc. 923 East Wellspring Road New Freedom, PA, USA 17349 +1 717-235-6807 www.oakworks.com	<b>Manuf. date:</b>	2023-06-27
	<b>SN</b>	CFA840463
	<b>MD</b>	QTY: 1 EA

(01)00817463022506(11)230627(21)CFA840463

This product complies with United States Department of Health and Human Services radiation performance standards, 21 CFR Subchapter J, in effect at the time of manufacture for Radiologic Patient Cradles.

**OAKWORKS, INC.**  
923 East Wellspring Rd.  
New Freedom, PA 17349 USA  
Phone: 717-235-6807  
www.oakworks.com

made in the USA  
with US & imported parts

**EC REP**

**EMERGO EUROPE**  
Westervoortbedijk 60  
6827 AT Arnhem, The Netherlands  
  
Part No. CN71047 - Rev. 2

### Carbon Fiber Arm Board-2in Pad

<b>Manufacturer</b>	<b>GTIN</b>	00817463022506
OAKWORKS, Inc. 923 East Wellspring Road New Freedom, PA, USA 17349 +1 717-235-6807 www.oakworks.com	<b>Manuf. date:</b>	2023-06-27
	<b>SN</b>	CFA840463
	<b>MD</b>	QTY: 1 EA

(01)00817463022506(11)230627(21)CFA840463

↑

GTIN

↑

Manufacture Date

↑

Serial Number

Model Name

Unique Device Information:

- GTIN - 14 digit number unique for each variation of a model
- Manufacture Date - Country of manufacture and date the device goes into production in YYYY-MM-DD format
- CE - CE Mark
- SN - Serial Number
- MD - Medical Device Symbol
- QTY - Quantity of the product

---

THIS PAGE IS INTENTIONALLY LEFT BLANK

# USER MANUAL

---

## OAKWORKS®

### Imaging and Pain Management Accessory - Carbon Fiber Armboard

#### CONTACT INFORMATION:

**OAKWORKS® Inc.**

923 East Wellspring Road

New Freedom, PA USA 17349

Phone: 717.235.6807

FAX: 717.235.6798

[www.oakworksmed.com](http://www.oakworksmed.com)

**CoYoMe B.V.**

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3261 LD Oud-Beijerland, The Netherlands


[www.coyome.eu](http://www.coyome.eu)

Enrico Cohen

[Enrico@coyome.nl](mailto:Enrico@coyome.nl)

[info@coyome.nl](mailto:info@coyome.nl)

Phone: +31 613.886.424

Emergo Europe  
Westervoortsedijk 60  
6827 AT Arnhem  
The Netherlands

Original Instructions  
Manual Part Number MMINMD0003-EN / Med-Eng-PM-84  
Revision Level: 8  
Revision Date: 04/01/2025

Printed in USA