

# OAKWORKS® WARRANTY

Effective February 16, 2009

TABLE	FRAME / UNDERSTRUCTURE	UPHOLSTERY / PADDING	ELECTRICAL & MOVING PARTS	REPLACEMENT SHIPPING	RETURN SHIPPING
<b>PROFESSIONAL LINE PORTABLES</b>					
Advanta, Alliance Aluminum, Alliance Wood, Aurora, BOSS, Nova, WellSpring II, Desktop Portal, PortalPro®, Nova with ABC system	Limited Lifetime-original owner.	5 yrs.	na ABC = 2 yrs.	included	included
	International 10 yrs.			International not included	International not included
<b>PORTABLES</b>					
One & KeLa™, Utopian with ABC system	Limited Lifetime-original owner.	2 yrs.	na ABC = 2 yrs.	not included	not included
	International 10 yrs.				
<b>LIFT TABLES</b>					
Celesta® Gemini, Clodagh Gemini, Celesta® Deluxe, Celesta® Lift, Celesta® Performa Lift, ProLuxe Veneto, ProLuxe Europa, ProLuxe Provence, ProLuxe Valencia, ProLuxe Seville, Articulating Top, Antigua, Antigua Colonial, Clinician Lift, Marina, Relaxation Lounge, Clodagh Leo, Clodagh Virgo, Clodagh Libra	Limited Lifetime-original owner.	5 yrs.	3 years parts 1 year labor	included	not included
	International 10 yrs.			International not included	
<b>MAXX™</b>					
MAXX™	4 yrs.	4 yrs.	3 years parts 1 year labor	included International not included	not included
<b>STATIONARY</b>					
Spa Clinician, Classic Clinician, Studio Clinician, Caban Lounger	Limited Lifetime-original owner.	5 yrs.	3 years parts 1 year labor	included	not included
	International 10 yrs.			International not included	
<b>PT/AT/ORTHOPEDIC</b>					
Storable Mat, Portable Taping Table, Portable Manipulation Table	Limited Lifetime-original owner.	5 yrs.	na	included	included
	International 10 yrs.			International not included	International not included
PT100, PT200, PT300, Casting Table	Limited Lifetime-original owner.	5 yrs.	3 years parts 1 year labor	included	not included
	International 10 yrs.			International not included	
Power Line, Seychelle, OAKWORKS® Prone Pillow	3 yrs.	3 yrs.	na	not included	not included
<b>MEDICAL</b>					
Fluoroscopy/C-Arm imaging tables, Ultrasound Table Line	4 yrs.	4 yrs.	3 years parts 1 year labor	included	not included
				International not included	
Spine Positioning System	2 yrs.	1 yr.	na	included	included
		International 2 yrs.		International not included	International not included
<b>ACCESSORIES</b>					
Arm Hammock, Adjustable Face Rest, Carry Cases, Basic Stools, Table Fleece, Face Rest Cozy, Sheets, Heating Pad, Hot Towel Cabi, OAKWORKS® Trolley, Table Cart	1 yr.	1 yr.	na	not included	not included
All Other Accessories	2 yrs.	2 yrs.	na	included	included
				International not included	International not included

Please register your table warranty online. [WWW.OAKWORKS.COM](http://WWW.OAKWORKS.COM)

# How to Care for Your Upholstery Fabric

Congratulations on your table purchase. In order to maintain your product's quality and ensure the life of your table, we have tested several preventative steps you can take to protect the upholstery fabric on your table and accessories.

Here are a few suggestions in caring for your upholstery fabric:

- DO NOT store your table or upholstery fabric-covered accessories in temperatures below 32°F/0°C or above 95°F/35°C. Extreme hot and cold temperatures will damage the upholstery fabric on your table and accessories. Some potential results are cold-cracking, brittleness and upholstery fabric stretching. Be sure to take your table out of the trunk of your car during extreme heat or cold. And, always keep your table away from heaters, electrical appliances and direct sunlight.
- Protect your table and accessories with sheets and fleeces. We offer both fitted sheets and thick fleeces that protect your upholstery fabric against oil stains while offering comfortable 100% cotton flannel softness to your clients. You can also preserve your face rest and bolsters with 100% cotton flannel covers and fleeces.
- Use a carrying case. We offer carrying cases to protect your table against scuffs, scrapes and punctures when travelling.
- Resist oil stains with the Anti-Slip Pad Our Anti-Slip Pad serves two functional purposes--it guards your table's upholstery fabric from oil that may damage the upholstery fabric over time, and it provides a non-slip surface under the sheet on your table. The Anti-Slip Pad is available in 24"x 65".
- Clean your upholstery fabric every day. We recommend cleaning your upholstery fabric with a mild detergent to protect it against oil stains and dirt which will cause cracking and permanent stains over time. Household products you can use include a 4:1 diluted solution of 409® cleaner, Fantastik® or any non-abrasive with a non-alcohol content. Products with a high alcohol content, when used over time, will dry and eventually crack your upholstery fabric.
- Use both ends of the table evenly. Occasionally switch the face rest from one end of the table to the other to reduce uneven wear.

Caring for your table's upholstery fabric involves more than just protecting your table from general abrasions, cuts and punctures. We offer most of the products mentioned above in our therapeutic equipment catalog.

Tables may be reupholstered for a charge, but you can lengthen the life of your upholstery fabric by following these quick and easy suggestions.

*If you have any questions, call Customer Service  
at 717-235-6807 or visit our website:  
[www.oakworks.com](http://www.oakworks.com)*

